


Generic Respirator Mask Sizing Chart

Note: *Approximately 75-80% of people fit into the medium size respirator*



We recommend taking exact measurements using the method described below to determine the correct mask size for your facial structure.

Measure from the centerline of your face at your mouth to at least 1 inch (2.5 cm) past the outside corner of your mouth (you might want to ask a friend to assist you). The outer edge of the Mask should sit against your face at 1 to 1.5 inches (2.5 to 3.8 cm) outside the corner of your mouth.

Double the measurement you have taken and match to the best fit below:

LARGE - 5 inches (12.7 cm)

MEDIUM - 4.5 inches (11.4 cm)

SMALL - 4 inches (10.2 cm)